

COVID-19 ECOTEST

ANTIGEN & ANTIBODY RAPID TEST





ANTIGEN RAPID TEST



NASOPHARYNGEAL SWAB

The COVID-19 antigen rapid test device is a lateral flow in vitro immunoassay intended to be used for the qualitative detection of the SARS-

CoV-2 viral nucleoprotein. This device will detect both viable and nonviable virus in both nasopharyngeal and oropharyngeal secretions.

ANTIBODY RAPID TEST




FINGERSTICK WITH LANCET

Detect the presence of SARS-CoV-2 IgM and IgG antibodies instantly from a whole blood sample collected from a minimally-invasive fingerstick blood collection process.

ANTIGEN & IgM/IgG DETECTABLE TEST WINDOW

The COVID-19 ECOTEST Antigen Rapid test will detect the antigen for 5 days after symptoms first appear.

The COVID-19 ECOTEST Antibody Rapid test will detect antibodies 5 to 7 days after symptoms first appear.







PROCEDURES AND READING RESULTS

ANTIGEN STEPS









Т





NEGATIVE: Only one colored

band appears, in the control

No apparent colored band

appears in the test region (T).

region (C).





Drop






Read

RESULT INTERPRETATION



C



One band appears in the control region (C) and another band appears in the test region (T).



INVALID: Control band fails to appear.

Results from any test which has not produced a control band at the specified read time must be discarded. Please review the procedure and repeat with a new test. If the problem persists, discontinue using the kit immediately and contact your local distributor.

NOTES:

- The colour intensity in the test region (T) may vary depending on the concentration of analytes present in the specimen. Therefore, determine the concentration of analytes in the specimen.
- Insufficient specimen volume incorrect operating procedure or expired tests are the most likely
- This test is intended for

Positive Percent Agreement (PPA): 94.3 % (95% Cl:84.6% - 98.1%) Negative Percent Agreement (NPA): 98.3 % (95% CI:95.1% ~ 99.4%)


Overall Agreement: 97.4 % (95% CI:94.4% ~ 98.8%)

ANTIBODY STEPS

C

т









Pipette blood sample





Add buffer, read results

IgG and IgM POSITIVE

The coloured line in the control


line region (C) changes from

blue to red, and two coloured

lines appear in test line regions


The result is anti-COVID-19 IgM

RESULT INTERPRETATION





С

IgG









NEGATIVE

The coloured line in the control line region (C) changes from blue to red. No line appears in the test line regions M or G.

The result is Negative.

IgM POSITIVE

The coloured line in the control line region (C) changes from blue to red, and a coloured line appears in test line region M.

The result is anti-COVID-19 IgM Positive.

IgG POSITIVE

The coloured line in the control line region (C) changes from blue to red, and a coloured line appears in test line region G.

The result is anti-COVID-19 IgG Positive.

INVALID

and IgG Positive.

Control line is still completely or partially blue, and fails to completely change from blue volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new cassette. If the problem persists, discontinue using the kit immediately and contact your distributor.

IgM: Relative sensitivity: 93.7% (86.0%-97.3%)*
Relative specificity: 99.1% (96.8%-99.8%)*
Overall agreement: 97.7% (95.4%-98.9%)*
*95% Confidence Interval

IgG: Relative sensitivity: 98.8% (93.5%-99.8%)*
Relative specificity: 98.7% (96.2%-99.5%)* Relative specificity: 98.7% (96.2%-99.5%)*

Overall agreement: 98.7% (96.7%-99.5%)*
*95% Confidence Interval



BLUE

RED

RED

ORDER YOURS TODAY.

asmith@promedica.ca

CELL: 780-913-5277 | OFFICE: 780-455-0368